

## 2.3 Software Solutions

### 2.3.1 StarLab

StarLab turns your PC into a laser power/energy multi-channel station

#### Extensive Graphic Display of Data

- Line Plot, Histogram, Bar chart, Simulated Analog Needle
- Multiple data sets on one graph or separate graphs on the same screen

#### Advanced Measurement Processing

- Power/Energy Density, Scale Factor, Normalize against a reference
- Multi-channel comparisons
- User defined mathematical equations: channels A/B, (A-B)/C etc.
- Position & size measurement with BeamTrack sensors

#### Data Logging for Future Review


- Can be displayed graphically or saved in text format
- Easily exported to an Excel spreadsheet

Fully supports Centauri, StarBright, StarLite, Vega, Nova-II, Pulsar, Juno, Juno+, Quasar, EA-1 and USB1 devices with all standard Ophir sensors


#### Flexible Display Options with StarLab

You may choose to display them separately


Choose which channels to display




Maximize one of the sources




Setup screen



Choose line graph




or histogram



One of the above screens is maximized


or needle display

Multiple Sensors displayed together



Here multi line graph display has been chosen

Settings and functions may be opened to adjust then minimized as needed



Here multi line histogram display has been chosen

## Functions and Logging

### Functions

Click on f(x) to open another trace combining measured values

Define function combining measured values

New trace is now added per defined function

### Logging

Click on log button and logging of values starts

Files are stored here. They may be viewed graphically OR numerically


```

:PC Software:StarLab version 3.00 Build 19
:Logged:25/05/2014 at 09:33:22
:Channel B:vega Thermopile 3A-P-V1 (s/n:999999) V2.31 (s/n:657028)
:Channel A:Juno Photodiode PD300 (s/n:694646) 3N1.24 (s/n:606180)
:Math M:(A-B)^2
:Channel B:Statistics
:Min:3.440mw
:Max:12.22mw
:Average:7.882mw
:Std.Dev.:3.078mw
:Overrange:0
:First Pulse Arrived : 25/05/2014 at 09:33:22.562000

```

Timestamp	Channel B	F(B)	Channel A	Math M
0.000	1.762e-002	6.620e-003		
0.064	1.836e-002	7.390e-003		
0.128	1.911e-002	8.110e-003		
0.136			1.067e-002	6.554e-006
0.193	1.986e-002	8.860e-003	8.480e-003	1.444e-007
0.203			6.540e-003	9.181e-006
0.256	2.057e-002	9.570e-003	4.900e-003	2.841e-005
0.269	2.123e-002	1.023e-002	3.550e-003	5.285e-005
0.321	2.182e-002	1.082e-002		
0.354			3.400e-004	1.339e-004
0.384	2.232e-002	1.132e-002	3.600e-004	1.259e-004
0.406	2.291e-002	1.191e-002	4.800e-004	1.141e-004
0.449	2.258e-002	1.158e-002	7.600e-004	9.761e-005
0.865	2.216e-002	1.116e-002	1.340e-003	7.569e-005
0.870			2.370e-003	4.914e-005
0.928	2.164e-002	1.064e-002		
0.936	2.104e-002	1.004e-002		
0.993	2.038e-002	9.380e-003		
1.003				
1.056	1.558e-002	4.580e-003		
1.070				
1.120				
1.136				
1.184				
1.203				
1.664				

# BeamTrack Power/Position/Size Screens



Power / Position / Size screen



Position stability screen

Displays beam center wander weighted for dwell time at each position