

# How to easily speed up end-of-line quality testing for lasers

Authors: Christian Dini and Sven Schipper

The laser is a very powerful and universal tool. It has found its way into many processes across all industries, including medical applications where the number of use cases continues to grow. Reliability and safety demands are placing high levels of responsibility on manufacturers of laser systems. To ensure the high quality of each laser system, thorough end-of-line testing is mandatory before a laser system is delivered. Nevertheless, as with many industries today, competition is high and manufacturers are seeking to more cost-effectively produce while meeting high-quality standards. Using the Ophir® BeamSquared® system, end-of-line-testing for laser systems can be performed much faster and the process can be optimized.

## *M2 and the quality of a laser*

When it comes to determining the quality of a laser, the key parameter is the beam propagation ratio or  $M^2$  (M-Squared). Its value indicates how close the laser beam comes to be a single mode TEM00 beam, which in turn determines how small a beam waist can be focused. For the perfect Gaussian TEM00 condition,  $M^2$  equals 1. While M-Squared is a key parameter for laser manufacturers, it usually requires some measurement efforts: To determine the beam propagation ratio, as described in the ISO/DIS 11146 standard, a series of measurements are mandatory.  $M^2$  is measured on real beams by focusing the beam with a fixed position lens of known focal length, and then measuring the characteristics of the artificially created beam waist and divergence.




Figure 1: Beam width measurements of a focused beam.

To provide an accurate calculation of  $M^2$ , it is essential to make at least five measurements in the focused beam waist region and at least five measurements in the far field, two Rayleigh ranges away from the waist area. The multiple measurements ensure that the minimum beam width is found. In addition, multiple measurements enable a "curve fit" that improves the accuracy of the calculation by minimizing measurement error at any single point. An accurate calculation of  $M^2$  is made by using the data from the multiple beam width measurements at known distances from a lens, coupled with the known characteristics of the focusing lens.

