

ModeCheck Readme  
Version 6.13

**Contents:**

- Section 1 - Change log
- Section 2 - Errata and Workarounds
- Section 3 - ModeCheck Notes

\*\*\*\*\*

**Section 1 - Change log**

V6.13.0 04/23/2018

- Upgraded FlyCap camera driver to 2.12.2
- Corrected an error in display of beam intensity color values in the Color Bar, 2D cursor tooltip and status bar when power calibrated.
- Improved performance of most results.

V6.12.0 02/26/2018

- Corrected an error in the calculation of Positional Stability.

V6.11.0 9/13/2017

- Upgraded PGR drivers to 2.11.164.
- Added support to auto upgrade installed drivers during installation.
- Improved reliability of the console service and data server interaction allowing connections to devices to be more stable.

V6.8 - 6/28/2016

- Full Windows 10 compatibility.
- Enhanced application logging for diagnostics.
- Renamed application title to make it easier to identify in Task Manager.

\*\*\*\*\*

**Section 2 - Errata and Workarounds**

We work hard to find and correct any bugs in this software product. However, as of this release we still have a few tough bugs for which we have not found complete solutions. The following list details these bugs and offers recovery and work-around methods if available:

- Users with Adobe Reader X may see inconsistent behavior when using the What's This? feature. This is due to flaws in a new security feature in Adobe Reader X. Users can restore previous functionality of What's This? by disabling "Protected Mode" in the Reader X Preferences menu.

\*\*\*\*\*

**Section 3 - ModeCheck Notes**

Supported Operating Systems:

- Windows 10 (64 bit)
- Windows 7 (32-bit & 64-bit)
- Windows Vista (32-bit & 64-bit)

While ModeCheck is a Windows Vista/7/10 product, full ModeCheck operation is expected in earlier Windows operating systems. Extensive testing has not and will not be performed in earlier operating systems; however, we have yet to encounter any major operational problems.

**Documentation:**

- A PDF version of the Operator's Manual is included with the installation. You must have Adobe Acrobat Reader in order to view this file. You can install Adobe Acrobat Reader from <https://get.adobe.com/reader/>.

**Installation:**

- It is recommended that all users are fully updated to the latest Windows Updates. If all updates are not applied to your system this may cause problems with ModeCheck.
- Under Windows Vista/7/10 you must have Administrative privileges in order to fully install ModeCheck and the required camera driver package.

If you suspect you have found a bug in our software please help us identify it by sending a description of the actions that reproduce it and the .mcSetup file you were using at the time to [softwareqa@ophir-spiricon.com](mailto:softwareqa@ophir-spiricon.com). The more information you can provide, the more likely we can reproduce it in our lab, and fix it.

*\* ModeCheck is a registered trademark of Ophir-Spiricon, LLC*

*\* Windows, Windows XP, Windows Vista, Windows 7, and Windows 10 are registered trademarks of Microsoft Corporation in the United States and other countries.*