

Ophir-Spiricon Camera Types That Support Power/Energy Interface With StarLab:										(Rev 12)	12/19/2011				
Cameras	SP620	GRAS20	L11058	LW230	Xeva	ModeCheck	BeamCube	Pyrocam III	Pulsar	Vega	Nova II	Quasar	Juno	LaserStar	
	SP503	SCOR20							USBI					Single/Dual	
	(USB)	(FW)	(USB)	(USB)	(USB)	(USB)	(USB)	(FW)	(1,2,&4)					Nova/Orion AN2	
BeamGage 32 bit Windows OS	OK	OK	OK	OK	OK	OK	OK	OK							
BeamGage 64 bit Windows OS	OK	NO	OK	OK	NO	NO	OK	OK							
StarLab 32 & 64 Windows OS									OK	OK	OK	OK	OK	NO	
BeamGage and StarLab 32 bit OS	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	NO	NO	NO	
BeamGage and StarLab 64 bit OS	NO	NO	NO	NO	NO	NO	NO	NO	NO	NO	NO	NO	NO	NO	
BeamGage Compatible		StarLab Compatib		Not Compatible				Not Recommended							

Photon 64bit vs 32bit Compatability Chart:

(Rev 12)

12/19/2011

Photon Equipment	Nano Scan <i>(USB)</i>	Nano Scan <i>(PCI)</i>	Nano Mode Scan <i>(USB)</i>	Nano Mode Scan <i>(PCI)</i>	M ² 1780 <i>(FW)</i>	Gonio LD8900 <i>(8 Bit)</i>	GonioR LD8900R <i>(16 Bit)</i>
------------------	---------------------------	---------------------------	--------------------------------	--------------------------------	------------------------------------	--------------------------------	-----------------------------------

Windows OS 32 bit	OK	OK	OK	OK	Vista OS or Earlier	Vista OS or Earlier	Vista OS or Earlier
-------------------	----	----	----	----	---------------------------	---------------------------	---------------------------

Windows 64 bit OS	OK	NO	OK	NO	NO	NO	NO
-------------------	----	----	----	----	----	----	----

Compatible

Limited Compatability

Not Compatible