












Display Finder

The table below lists the specs and features of Ophir Smart Displays and PC Interfaces

	Displays						PC Interfaces				Wireless Interface
	<u>Orion</u>	<u>Nova</u>	<u>Nova II</u>	<u>Vega</u>	<u>LaserStar Single Channel</u>	<u>LaserStar Dual Channel</u>	<u>USB1</u>	<u>Pulsar-1</u>	<u>Pulsar-2</u>	<u>Pulsar-4</u>	<u>Quasar</u>
											
Digital Display	Yes	Yes	Yes	Yes	Yes	Yes	N/A	N/A	N/A	N/A	N/A
Display Color	Monochrome	Monochrome	Monochrome	Color	Monochrome	Monochrome	N/A	N/A	N/A	N/A	N/A
Analog Display	No	No	Yes	Yes	No	No	N/A	N/A	N/A	N/A	N/A
Rechargeable Battery	Yes	Yes	Yes	Yes	Yes	Yes	Bus Powered	No	No	No	Yes
Detector Support:											
Thermal Heads	Orion-TH	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Photodiode Heads	Orion-PD	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Pyroelectric Heads	Orion-PE	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
RP Heads	No	No	No	No	Yes	Yes	Yes	No	No	No	No

