

## Main Memory Volatility Statement

---

### **Models:** Photodiode Sensors

#### **Product Description:**

General Purpose Laser Power Sensors

#### **Memory Description:**

These sensors contain the following memory devices:

##### Smart plug:

U3: I2C EEPROM 93C66,

Used to store sensor parameters and calibration factors. Non-volatile.

#### **General:**

Sensor calibration constants are stored in U3 in smart plug board. The calibration constants are generated when the sensor is sent through its calibration process in the factory, and are fundamental to the sensor operation.

**Note:** The standard sensors contain a D15 plug for power meter connection. The measuring of the sensors displayed on a standard Ophir meter or PC interface.

Gadi Seeberger  
QA manager